

OxiMax™ Single Patient Use Sensors



Type of Patient	Adult	Paediatric	Infant	Neonate	Adult	Premature Neonates	Neonate	Adult	Adult/ Paediatric
Weight Range	> 30 kg	10 to 50 kg	3 to 20 kg	< 3 kg	> 50 kg	< 1.5 kg	1.5 to 5 kg	> 40 kg	> 10 kg
Sensor Type	MAX-A-I MAX-AL-I	MAX-P-I	MAX-I-I	MAX-N-I	MAX-R-I	Softcare™ SC-PR-I	Softcare™ SC-NEO-I	Softcare™ SC-A-I	MAX-FAST™-I
SpO ₂ Range LoSat™	70% - 100%	± 2	± 2	± 2	± 2	± 3.5	± 2	± 2	± 2
	60%-80%	± 3	± 3	± 3	± 3				± 3

OxiMax™ Reusable Sensors & Accessories - Wraps



Type of Patient	All	Adult	Adult/ Neonate	Paediatric/ Infant	Adult/ Paediatric	Type of Patient	All	Adult/ Neonate	Adult/ Neonate	Paediatric/ Infant	Paediatric/ Infant	
Weight Range	> 1 kg	> 40 kg	< 3 kg or > 40 kg	3 to 40 kg	> 30 kg	Weight Range	All	< 3 kg or > 40 kg	< 3 kg or > 40 kg	3 to 40 kg	3 to 40 kg	
Sensor Type	D-YS	DS100A	OXI-A/N	OXI-P/I	D-YS with D-YSE	Sensor Type	D-YS	OXI-A/N	OXI-A/N	OXI-P/I	OXI-P/I	
SpO ₂ Range LoSat™	70% - 100%	± 3*	± 3	± 3*	± 3	± 3.5	P/n fixation	Posey	FOAMA/N	ADH-A/N	FOAMP/I	ADH-P/I

* Accuracy in neonates is ± 4

Accuracy ranges stated have been validated with the N600x OxiMax™ monitor